All Purpose Durable BP Cuff

English

Suntech Medical's All Purpose Clinical Grade Blood Pressure Cuffs

INDEX  Index Line

 Cuff index line must fall within range markings

ARTERY  Artery symbol and arrow should be placed over brachial or femoral artery

 Not made with PVC

 Not made with natural rubber latex

CE  Product in compliance with Council Directive 93/42/EEC Medical Devices Directive

 Symbol indicating arm circumference

 Symbol indicating manufacturer

LOT  Symbol indicating lot code of cuff

 Symbol indicating Caution

RANGES / COLORS

Size:	Color:	Range:
Infant	Orange	8-13 cm
Child	Green	12-19 cm
Child Long	Green	12-19 cm
Sm Adult	Royal Blue	17-25 cm
Sm Adult Long	Royal Blue	17-25 cm
Adult	Navy Blue	23-33 cm
Adult Long	Navy Blue	23-33 cm
Lg Adult	Burgundy	31-40 cm
Lg Adult Long	Burgundy	31-40 cm
Thigh	Brown	38-50 cm

WARNINGS

Pressurization of the cuff can temporarily cause loss of function if simultaneously used with monitoring equipment on the same limb.

The cuff should not be placed on the patient’s arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.

If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

Do not attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.

Following the application of the BP cuff, petechia formation (a minute reddish or purplish spot containing blood that appears in the skin) or Rumpel-Leede phenomenon (multiple petechia) on the arm, which may lead to idiopathic-thrombocytopenia (spontaneous persistent decrease in the number of platelets, associated with hemorrhagic conditions) or phlebitis (inflammation of a vein) may be observed.

CAUTIONS

After washing ensure the size indication on the bladder and cuff shell match. Make sure that the cuff hose is threaded through one of the hose openings in the cuff.

Minimize limb movement during the measurement.

The cuff should not be applied over a wound as this can cause further injury

Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner
Do not use cuff if there are any signs of damage. Failure to do so could affect measurement accuracy

Do not machine wash the cuff bladder. Water could be trapped in the cuff and cause damage to the NIBP module and/or inaccurate BP readings
Avoid contact with the cuff, other than that of the patient’s limb, while measurement is in progress. A compressed or kinked connection hose may

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SunTech Medical, Inc.
SunTech Medical, Ltd. Europe