



STAINLESS



MODEL 7748SS (STANFORD) MAYO STAND maneuvers easily with two additional casters attached to a "U" shaped tubular frame. Great for the ER or exam rooms. Height adjusts with a Tru-Loc™ friction knob. Also available with optional MODEL 7748CR CRADLE ASSEMBLY, including a carriage to house a 13 qt. kickbucket (kickbucket included).

Model No.	Part No.	Feature	Height Adjustment	Tray Size	Ship Wt. (lb) est.
7748SS	0647748000	"U" Shape base	39.5 "- 62"	13-3/4" x 19-1/8"	35

Option: 12 qt. kickbucket inserted into frame.

Option: 7748CR CRADLE ASSEMBLY, includes carriage to house 13 qt. kickbucket (kickbucket included).







MODELS 7767SS AND 7769SS (RICHMOND) FOOT OPERATED MAYO STANDS are number one in ease and convenience. Just lift the tray frame – it automatically locks. To lower, hold the tray frame and press the foot operated control button. Release to lock in place. A guard over the foot button prevents accidental release. Select our MODEL 8771SS, DOUBLE-POST JUMBO MAYO STAND when a larger tray size is required. MODELS 7740SS AND 7749SS (MANHATTAN) MAYO STANDS lower with a thumb-operated control button. To raise, lift the tray frame to desired height – it automatically locks in place.

Model No.	Part No.	Feature	Height Adjustment	Tray Size	Ship Wt. (lb) est.
7767SS	0617767000	Foot-operated	39.5" - 62"	13-3/4" x 19-1/8"	46
7769SS	0617769000	Foot-operated	39.5" - 62"	16" x 21"	46
8771SS	0658771000	Foot-operated	39.5" - 62"	20" x 25"	54
7740SS	0627740000	Thumb-op. button	39.5" - 62"	13-3/4" x 19-1/8"	45
7749SS	0627749000	Thumb-op. button	39.5" - 62"	16" x 21"	46

